

 **cuattro DR^{HD}**




exclusive



A whole new category of Digital Radiography...

Cuatro DR HD (High Definition) Digital Radiography is now available exclusively from Heska. Cuattro DR HD delivers detail, resolution, and contrast image quality that is visibly better, much better, than Standard Definition systems.

 **cuattro DR^{HD}**



HD Minimum Standards:

- Direct Deposition Cesium
- 100 micron pixel size
- HD pixel density. 18 Million (17x17)
- HD pixel density. 15 Million (14x17)
- Full Size Image Area (14x17 - 17x17)
- Resolution of at least 5 lp/mm
- 16 Bit A/D 16 Bit Grayscale Native
- 16 Bit Multi-Frequency Tuning



Cuattro DR HD Added Benefits:

- Automatic Exposure Detection
- Cassette Sized Form Factor
- Embedded - 64 bit OS
- Multi-Touch User Interface
- Cloud Based PACS Connectivity
- Cloud Archival
- Software Upgrades and Updates
- 24/7 Live Technical Support

A smarter direction in Digital Radiography. Better, much better.

www.heska.com • toll free: 800.709.4515 • imaging the possibilities

© 2019 Heska Corp. All rights reserved. Specifications are subject to change without notice and may vary based on operating conditions. Names of companies or products appearing in this document are trademarks of their respective owners. Cuattro DR is a trademark of Heska Corporation in the United States and may also be a registered trademark in other countries. 19LT0105

