



Introducing the Optomed VetCase™ portable video endoscope.

Designed with veterinarians in mind, the Optomed VetCase system integrates the insufflation and processor unit, LED light source, and video and acquisition components into a portable and rugged self-contained package. The large ergonomic touchscreen interface displays and records video and stills for real-time viewing.

The Optomed VetCase accepts all types of rigid and video endoscopes and is perfect for both equine and small animal practices. Other Optomed scopes can be used, including the fiberscope, otoscope, arthroscope and laparoscope. This flexibility allows for a variety of applications, such as equine larynx, guttural pouches, insemination, bronchi, and gastroscopy. Full studies can be integrated into CloudBank for anywhere, anytime, any device viewing with Web access.

www.heska.com • toll free: 800.709.4515 • imaging the possibilities

© 2019 Heska Corporation. All rights reserved. HESKA is a registered trademark and Optomed VetCase and CloudBank are trademarks of Heska Corporation in the United States and in other countries. Specifications are subject to change without notice and may vary based on operating conditions. Pricing dependent on final configuration. US19LT0608