



Now with  
Progesterone,  
cTSH & cCRP

# element i+

Immunodiagnostic Analyzer

## The Future of Immunodiagnostics is Here!

Combine proprietary technology with a sleek, compact platform to bring you fast, accurate results with an expanding test menu.

Fast results in as little as 4 minutes

Room temperature cartridge storage

Accurate, quantitative results

Large color touch screen

Bi-directional interface

## Easy Operation 3-step easy measurement



1 Enter Patient Information



2 Load Patient Sample



3 Insert Cartridge

### Features:



## Test Parameters

### MEASURED PARAMETERS

|              |  |
|--------------|--|
| T4           | Diagnose and monitor thyroid disorders                               |
| cTSH         | Diagnose and monitor hypothyroidism in canines                       |
| Cortisol     | Diagnose and monitor patients with Cushing's or Addison's disease    |
| Progesterone | Manage breeding and parturition and screen for reproductive diseases |
| cCRP         | Screen and monitor systemic inflammation in canines                  |

### TESTS IN THE PIPELINE

Cancer screening, infectious disease, SDMA, and more.

### Convenience

- Fully self-contained micro-fluidic cartridge
- Simplified interface and operation with a small bench top footprint
- Small 100µL sample volume of serum or plasma

### Accuracy

- Fluorescence immunoassay provides accurate, quantitative results
- Superior accuracy of patented, proprietary fluorescence waveguide immunoassay technology with laser evanescent illumination

## In-Clinic

### Expand your in-house diagnostic capabilities

- With run times of 10 minutes and under, you can diagnose patients and develop treatment plans within the first visit
- No more waiting for lab results and playing phone-tag with clients
- Quickly and easily schedule follow-up monitoring with complex patients

## Results and Reporting

### Efficient data management

- Integrate with Practice Management Software
- Bi-directional communication enabled auto patient entry and medical record reporting
- Sort options make retrieving previous patient results a snap
- Create comprehensive patient reports with results from other Heska analyzers



www.heska.com 800.464.3752

©2022 Heska Corporation. All Rights Reserved. HESKA and Element i+ are registered trademarks of Heska Corporation in the U.S. and other countries. US19LT1001\_005

