

NEW



# element i+

Immunodiagnostic Analyzer

## The Future of Immunodiagnosics is Here!

Combine proprietary technology with a sleek, compact platform to bring you fast, accurate results with an expanding test menu.

Fast results in as little as 5 minutes

Room temperature cartridge storage

Accurate, quantitative results

Large color touch screen

Bi-directional interface

## Easy Operation 3-step easy measurement



1 Enter Patient Information



2 Load Patient Sample



3 Insert Cartridge

### Features:

# Ei+

## Test Parameters

### MEASURED PARAMETERS

T4	Diagnose and monitor thyroid disorders
Cortisol	Diagnose and monitor patients with Cushing's or Addison's disease

### COMING SOON

Progesterone, TSH (Canine), Infectious Disease, and more

### Convenience

- Fully self-contained micro-fluidic cartridge
- Simplified interface and operation with a small bench top footprint
- Small 100 $\mu$ L sample volume of serum or lithium heparin plasma

### Accuracy

- Fluorescence immunoassay provides accurate, quantitative results
- Superior accuracy of patented, proprietary fluorescence waveguide immunoassay technology with laser evanescent illumination

## In-Clinic

### Expand your in-house diagnostic capabilities

- 5 minute run-time means you can diagnose and develop treatment plans on the first visit
- No more waiting for lab results and playing phone-tag with clients
- Quickly and easily schedule follow-up monitoring with complex patients

## Results and Reporting

### Efficient data management

- Integrate with Practice Management Software
- Bi-directional communication enabled auto patient entry and medical record reporting
- Sort options make retrieving previous patient results a snap
- Create comprehensive patient reports with results from other lab analyzers (Chemistry, Hematology, and others)



[www.heska.com](http://www.heska.com) 800.464.3752

©2020 Heska Corporation. All Rights Reserved. HESKA and HeskaView Integrated Software are registered trademarks and Element i+ is a trademark of Heska Corporation in the U.S. and other countries. US19LT1001

