
Article Information Sheet

Producer: FUJIFILM Corporation
26-30 Nishiazabu 2-Chome, Minato-ku, Tokyo 106-8620, Japan
Distribution: FUJIFILM Medical Systems U.S.A., Inc.; 419 West Avenue, Stamford, CT 06902;
Tel.: (203) 3242000, Fax: (203) 3530926;
Issued: 5/25/2018
Modified:

1. Product

Product Name: FUJI DRI-CHEM IMMUNO AU CARTRIDGE v-PRG

2. Composition

The main ingredients are:

Main Ingredients:	Weight-%:
acryl resin	70 - 80
polypropylene	20 - 30

3. Health & Safety Hazards

If used as intended, this product does not have any special health or safety concerns.
Take precaution for the bodily fluid samples being tested.

4. Fire Hazards and Fire Fighting Measures

This product is capable of burning, but will not ignite or flash under normal conditions.

Product / Main Ingredient(s):	Ignition Temp.	Flash Point	Density
product	n/av	n/av	n/av
acryl resin	450-462 °C	280-300 °C	n/av
polypropylene	>400 °C	340-360 °C	n/av

Combustion Gases:

Carbon dioxides, carbon monoxides, small amounts of nitrogen oxides, and smoke can be generated depending upon the burning condition.

Fire Fighting Measures:

Water, carbon dioxide, extinguishing powder or foam can be used as an extinguishing media.

While fire fighting, wear protective equipment such as self contained breathing apparatus, depending upon the fire situation.

5. Disposal Considerations

5.1 General Disposal Considerations

It is recommended to provide the waste to a legitimate recycler. If this is not possible, dispose the material in accordance with federal, state and local solid waste disposal laws.

5.2 Special Disposal Information for the US

After use, this product is categorized as an infectious waste. Dispose the used product in accordance with waste disposal regulations and use proper means of disposal, such as, incineration, melting, sterilization or disinfection.

6. Transport

This product is not subject to dangerous material classification and labelling pertaining to transportation and storage. The used product may be classified as a regulated Medical Waste, N.O.S. UN3291, if it has not been sterilized or otherwise treated.

7. Storage

Specific hazards regarding storage are not known. The product must be stored at 35.6-46.4 °F (or 2-8 °C) before use.

8. Further Information

Please refer to the operation manual for handling of this product.

The content of this information sheet has been drawn up in accordance with current knowledge at the time when the sheet was issued. This Article Information Sheet provides information related to the environment and to safety and does not make any warranty about the properties or quality of the product. Furthermore, the precautions presented herein are postulated upon normal handling (the type of use originally intended by the manufacturer) of the product.