

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0324A



2024-05-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0324A-1	LOT	ME0324A-2	LOT	ME0324A-3
element HT5+	Ret# x10 ³ /μL	88.2	± 64.4	229.5	± 99.0	181.7	± 99.0
QC Mode	Ret%	1.52	± 1.12	5.01	± 2.00	8.15	± 3.50
	IRF %	6.8	± 15.0	9.2	± 13.0	7.7	± 25.0
	LFR %	93.2	± 15.0	90.8	± 15.0	92.3	± 25.0
	MFR %	3.7	± 12.0	6.6	± 10.0	4.5	± 22.0
	HFR %	3.1	± 10.0	2.6	± 10.0	3.2	± 10.0
	RHE pg	28.2	± 12.0	24.6	± 12.0	23.1	± 12.0
	IPF %	7.7	± 15.0	7.6	± 15.0	7.6	± 15.0

Before using, refer to the instruction sheet for mixing directions.