

BC-RET

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT ME0923A



2023-11-10

Instrument	Parameter	Low		Normal		High	
		LOT	ME0923A-1	LOT	ME0923A-2	LOT	ME0923A-3
element HT5+	Ret# x10 ³ /μL	85.3	± 64.4	211.1	± 99.0	185.8	± 99.0
QC Mode	Ret%	1.47	± 1.12	4.67	± 2.00	8.22	± 3.50
	IRF %	6.4	± 15.0	8.2	± 13.0	8.6	± 25.0
	LFR %	93.6	± 15.0	91.8	± 15.0	91.4	± 25.0
	MFR %	2.7	± 12.0	5.7	± 10.0	6.2	± 22.0
	HFR %	3.7	± 3.7	2.5	± 3.0	2.4	± 3.0
	IPF %	6.6	± 15.0	5.2	± 15.0	5.5	± 15.0
	RHE pg	29.6	± 12.0	29.5	± 12.0	23.2	± 12.0

Before using, refer to the instruction sheet for mixing directions.